

Tuner-Gal

Product Name: Tuner-Gal

User Manual

Contents

Description	2
Application	2
Composition	2
Storage	2
General culture recommendations	2
Related Product	3

Description

Tuner-Gal is a glycosylation adjust with a clear chemical composition, which is independently developed, researched, and produced by Shanghai BioEngine Sci-Tech Co., Ltd. Tuner-Gal is designed to increase the galactosylation level of product protein in various cell cultures, such as Chinese Hamster Ovary (CHO) cells.

Application

This product is a sterile liquid product intended for research or further manufacturing in the biomanufacturing industry, but not for human or therapeutic use.

Composition

The IP rights of Tuner-Gal formulation are owned by Shanghai BioEngine Sci-Tech Co., Ltd.

This medium does not contain:

- Cytokines and antibiotics.
- Raw materials from animal sources.

Storage

- Store medium at 2-8°C, away from light.
- Do not use it after the expiration date or precipitation.

General culture recommendations

Tuner-Gal is added into the culture system at 0.1-0.2% (v/v, Tuner-Gal/Culture volume) on day 2 or 3 and then every other day.

Or

Add Tuner-Gal into basal medium with the volume ratio of 1:100 to 1:50 (*v/v*, Tuner-Gal/Basal medium volume).

Note:

- (1) The volume ratio could be further optimized according to the target glycosylation level.
- (2) Tuner-Gal can be used together with other Tuners to regulate protein glycosylation.

Related Product

Product Name	Туре	Cat. No.	Form	Size	Packaging	Note
Tuner-Gal	Additive	EXP0119201	Liquid	1 L	Bottle	Increase the galactosylation level of product protein in
		EXP0119202	Liquid	100 ml	Bottle	various cell cultures



Scan the QR code for more details about BioEngine.

Stay tuned for more updates.

Tel: 86-21-68582660

Web: www.bio-engine.com

