

Hyber-F100S Hybridoma Cell Feed Medium

Product Name: Hyber-F100S

User Manual

Contents

Description2
Application2
Composition2
Storage 2
Reconstitution of Powder Medium2
Fed-batch culture3
Related Product错误! 未定义书签。

Description

Hyber-F100S is a protein-free and animal-free feed medium specifically designed for fed-batch process of hybridoma cells. The medium enables the excellent growth performance of hybridoma cells and the high-level expression of antibodies, in conjunction with Hyber-B100 or Hyber-B100S medium (Refer to the "Related Product" section).

Application

Hyber-F100S feed medium is suitable for fed-batch culture with hybridoma cells.

This product is intended for research or further manufacturing in the bio-manufacturing industry, but not for human or therapeutic use.

Composition

The IP rights of Hyber-F100S feed medium formulation are owned by Shanghai BioEngine Sci-Tech Co., Ltd.

This medium contains:

- ✓ Carbohydrates, amino acids, vitamins, bulk salts, and trace elements.
- ☑ 28 g/L D-glucose, 2 g/L P188, 40 mM glutamine.

Not contain:

- ☑ Cytokines, antibiotics, HEPES and phenol red.
- \boxtimes Raw materials from animal sources.

Storage

- Store medium at 2-8°C, away from light.
- Once opened, the powder medium should be stored protected from moisture in a tightly sealed container.
- Do not use it after the expiration date or being damped.

Reconstitution of Powder Medium

Table 1 shows the preparation of Hyber-F100S medium [1].

Ingredients	Concentration
Hyber-F100S medium powder	123.18 g/L ^[2]
Sodium hydroxide	1.90 g/L
L-tyrosine	1.00 g/L
Hyber-B100S Additive	1 mL/L

Table 1. Preparation of Hyber-F100S medium

- Weigh 20% water of the final volume into the preparation container using pure water, ultrapure water, or water for injection at 20-30°C. Mix thoroughly (Power per Volume (P/V)>10 W/m³) without creating air bubbles.
- Accurately weigh 1.90 g/L of sodium hydroxide and add it into the preparation container of 1) step. Stir for 5-10 minutes.
- Weigh 1.00 g/L of L-tyrosine and add it into the preparation container of 1) step. Stir well for 5-10 minutes.
- Weigh the rest 80% water of the final volume into the preparation container of 1) step using pure water, ultrapure water, or water for injection. Accurately weigh 123.18 g/L of Hyber-F100S

medium powder into the preparation container. Stir well for 20-30 minutes.

- Adjust to pH 5.4-5.5 with 5 mL/L (recommended) concentrated hydrochloric acid solution. Stir for 20-25 minutes. At this point, the solution should be clear.
- 6) Add the Hyber-B100S Additive near the liquid level in the container.
- Slowly adjust to pH 6.0-6.3 with 5-10 mol/L sodium hydroxide solution. Stir for 10-15 minutes.
- Pass the medium solution through a pore size of
 0.22 or 0.2 µm sterile filter membrane, such as PES, using a pulse pump or compressed air (3-15 psi).
- 9) Use the prepared medium liquid immediately or store it in glass bottles, PET storage bottles, or disposable storage bags with an oxygen barrier membrane in a dark environment of 2~8°C. It's recommended for use within one month.

Note:

^[1] The above parameters (such as stirring time) are set for small-scale liquid preparation. Adjust these parameters for large-scale preparation based on container capacity to ensure full dissolution of dry powder.

^[2] The "g/L" unit denotes volumetric concentration (solute mass/solution volume).

Specifications of final liquid medium

Test	Unit	Specification
рН		6.0 – 6.3 ^[3]
Osmolality	mOsm/kg	600 – 1000
Turbidity	NTU	< 4.00

Table 2. Specifications of final liquid medium

Fed-batch culture

- The hybridoma cells should be adapted in Hyber-B100 or Hyber-B100S at least two passages.
- Ensure that the cell viability is >90%, and the growth rate is in mid-logarithmic phase prior to inoculation.
- Calculate the volume of cell culture and medium necessary to seed at 1.0-1.5×10⁶ viable cells/mL in a shake flask.
- Incubate at 37°C in a humidified atmosphere of 5% CO₂ in air on an orbital shaker platform rotating at 110-130 rpm (110 rpm for 50 mm amplitude; 130 rpm for 10 mm amplitude).
- The day of inoculation was recorded as day 0 (D0).
 Follow the suggested feed strategy as outlined in Table 3.
- Harvest the cells on day 7 or when viability falls below 50%.

Feed Medium	D1	D2	D3	D4	D5	D6	D7
Hyber-F100S Feed medium							
(%, v/v)	3	3	3	3	3	3	1

Table 3. Recommended feed strategy

Product	Cat. No.	Form	Size	Packaging	Notes
Hyber-B100 Hybridoma Cell Serum-free Medium	EXP0118901	Liquid	1 L	Bottle	• SF, PF, ADCF
Hyber-B100S Hybridoma Cell	EXP0111201	Powder	100 L	Bag	 Suitable for high-density
Serum-free Medium	EXP0111202	Powder	10 L	Bag	hybridoma cell suspension
	EXP0111203	Powder	5 L	Bag	cultures and supports high
Hyber-F100S Hybridoma Cell	EXP0111301	Powder	20 L	Bag	protein expression
Serum-free Feed Medium	EXP0111302	Powder	2 L	Bag	

Related Product



Scan the QR code for more product information.

Stay tuned for more updates. Tel: 86-21-68582660 Web: <u>www.bio-engine.com.cn/EN</u>

